

По вопросам продаж и поддержки обращайтесь:

Архангельск (8182)63-90-72	Краснодар (861)203-40-90	Рязань (4912)46-61-64
Астана (7172)727-132	Красноярск (391)204-63-61	Самара (846)206-03-16
Белгород (4722)40-23-64	Курск (4712)77-13-04	Санкт-Петербург (812)309-46-40
Брянск (4832)59-03-52	Липецк (4742)52-20-81	Саратов (845)249-38-78
Владивосток (423)249-28-31	Магнитогорск (3519)55-03-13	Смоленск (4812)29-41-54
Волгоград (844)278-03-48	Москва (495)268-04-70	Сочи (862)225-72-31
Вологда (8172)26-41-59	Мурманск (8152)59-64-93	Ставрополь (8652)20-65-13
Воронеж (473)204-51-73	Набережные Челны (8552)20-53-41	Тверь (4822)63-31-35
Екатеринбург (343)384-55-89	Нижний Новгород (831)429-08-12	Томск (3822)98-41-53
Иваново (4932)77-34-06	Новокузнецк (3843)20-46-81	Тула (4872)74-02-29
Ижевск (3412)26-03-58	Новосибирск (383)227-86-73	Тюмень (3452)66-21-18
Казань (843)206-01-48	Орел (4862)44-53-42	Ульяновск (8422)24-23-59
Калининград (4012)72-03-81	Оренбург (3532)37-68-04	Уфа (347)229-48-12
Калуга (4842)92-23-67	Пенза (8412)22-31-16	Челябинск (351)202-03-61
Кемерово (3842)65-04-62	Пермь (342)205-81-47	Череповец (8202)49-02-64
Киров (8332)68-02-04	Ростов-на-Дону (863)308-18-15	Ярославль (4852)69-52-93

Единый адрес: hno@nt-rt.ru **Веб-сайт:** www.hamilton.nt-rt.ru

**Емкости для хранения пробирок RackWare
Hamilton**

RackWare™

RackWare is Hamilton's new consumable product line featuring two types of SBS-footprint racks—standard racks and high-density racks. These racks are specifically designed to increase the storage capacity in our automated storage systems, but can also be used in manual freezers.

RackWare 48

This 48-position cryovial/microtube rack is a honeycomb design for use as an automation-friendly rack for larger labware. The rack can hold labware from Corning, FluidX, Greiner and Nunc™.



RackWare 48

RackWare HD60 and HD138

These high-density SBS-footprint racks for microtube and cryovial storage are designed for use with Verso, BiOS, and SAM systems. These racks increase the storage capacity in an automated or manual storage system by increasing the number of tubes that can be placed in the rack.



RackWare HD60

RackWare HD60 supports popular 48-format labware and increases storage capacity in a comparable SBS-footprint from 48 to 60 tubes.

RackWare HD138 supports most major 96-format labware and increases storage capacity in a comparable SBS-footprint from 96 to 138 tubes.



RackWare HD138

Technical Specifications:

	RackWare 48	RackWare HD60	RackWare HD138
Rack Dimensions (SBS Standard Footprint)	127.76 mm x 85.48 mm (5.03 in x 3.37 in)		
Temperature Range	-100°C to 120°C	-80°C to +35°C (-150°C in manual storage only)	-80°C to +35°C (-150°C in manual storage only)
Cryovial/Tube Diameter	10.0 mm to 12.0 mm (0.39 in to 0.47 in)	10.5 mm* (0.41 in)	8.5 to 9.0 mm* (0.34 in to 0.35 in)
Autoclavable	Yes	No	No

*Contact Hamilton for complete line of compatible tubes.



По вопросам продаж и поддержки обращайтесь:

Архангельск (8182)63-90-72	Краснодар (861)203-40-90	Рязань (4912)46-61-64
Астана (7172)727-132	Красноярск (391)204-63-61	Самара (846)206-03-16
Белгород (4722)40-23-64	Курск (4712)77-13-04	Санкт-Петербург (812)309-46-40
Брянск (4832)59-03-52	Липецк (4742)52-20-81	Саратов (845)249-38-78
Владивосток (423)249-28-31	Магнитогорск (3519)55-03-13	Смоленск (4812)29-41-54
Волгоград (844)278-03-48	Москва (495)268-04-70	Сочи (862)225-72-31
Вологда (8172)26-41-59	Мурманск (8152)59-64-93	Ставрополь (8652)20-65-13
Воронеж (473)204-51-73	Набережные Челны (8552)20-53-41	Тверь (4822)63-31-35
Екатеринбург (343)384-55-89	Нижний Новгород (831)429-08-12	Томск (3822)98-41-53
Иваново (4932)77-34-06	Новокузнецк (3843)20-46-81	Тула (4872)74-02-29
Ижевск (3412)26-03-58	Новосибирск (383)227-86-73	Тюмень (3452)66-21-18
Казань (843)206-01-48	Орел (4862)44-53-42	Ульяновск (8422)24-23-59
Калининград (4012)72-03-81	Оренбург (3532)37-68-04	Уфа (347)229-48-12
Калуга (4842)92-23-67	Пенза (8412)22-31-16	Челябинск (351)202-03-61
Кемерово (3842)65-04-62	Пермь (342)205-81-47	Череповец (8202)49-02-64
Киров (8332)68-02-04	Ростов-на-Дону (863)308-18-15	Ярославль (4852)69-52-93

Единый адрес: hno@nt-rt.ru **Веб-сайт:** www.hamilton.nt-rt.ru